Autoclave Indicator Tape

McKesson Autoclave Indicator Tape is used in securing wrapped instrument sets before, during, and after sterilization. Each roll of McKesson Autoclave Indicator Tape is individually packaged to ensure quality at the time of use. The Blue tape features a more aggressive adhesive and is for disposable wraps only.

Features and Benefits

- Adhesive for easy pack sealing
- Color change indicates that the pack was processed
- Color change is distinctly visible and non-reversible



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Autoclave Indicator Tape, Steam, Blue, 1” x 60 yds.	73-BT036	18ea/cs
Autoclave Indicator Tape, Ethylene Oxide (EO), White, 1” x 60 yds.	73-GT036	18ea/cs
Autoclave Indicator Tape, Steam, White, 1” x 60 yds.	73-ST036	18ea/cs
Autoclave Indicator Tape, Steam, White, ¾” x 60 yds.	73-ST048	24ea/cs

Satisfaction Guaranteed
If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care℠

With McKesson Brands, you’ll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates (“McKesson”). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson’s online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2015 McKesson Medical-Surgical Inc.

Rev. 00 08/15